

STREAM-4 FURNITURE SERIES



BONLAB
LABORATORY FURNITURE

STREAM 4

The STREAM-4 under bench furniture series is a robust furniture system on braked castors that can easily be moved around your laboratory

STREAM -4

Under bench units in this design are on castors and fully movable.

The configuration of furniture can be changed according to end user requirements by moving the under bench cabinets.

Mobile cabinets are fitted with worktop and can be used as separate bench by combining few cabinets together.

Front castors of each mobile units are with brake. Wheels are made of nylon wheels . Each wheel will carry 70kgs

The frame works of furniture benches have fine-adjust leveling footing.

**BONLAB STREAM-4 series
will provide quality and
style to any space.**

Our team's full experience has been applied to our benching systems – engineered to meet the demanding specifications of modern laboratories in all sectors.

