

STREAM - 2

Furniture structural series



STREAM- 2 series is designed with hygiene in action, Worktops are supported on Metal H' Frame (made of 30 / 40 mm square tubes) which provides easy cleaning of floors. This series is an ideal solution for hospitals or infection-control labs.

STREAM - 2


Furniture structural series



Al frames are welded to withstand high bearing loads, and also frames are epoxy powder coated for protection against chemical spills and corrosion.

Welding joints are easy cleanable and made with precision to avoid any possible contamination.

The frame works of furniture benches have fine-adjust leveling footing.



BONLAB **STREAM- 2** series system can be reinforced or adapted for heavier load bearing purposes

Our team's full experience has been applied to our benching systems – engineered to meet the demanding specifications of modern laboratories in all sectors.

