

BONLAB  
**STREAM**  
**1**  
SERIES

**STREAM-1 FURNITURE**  
SERIES



**STREAM-1 Series furniture range is a fixed system. Worktops will be fixed on cabinets and the Euro-style under bench cabinets are plinth mounted and does not require any framework**

# STREAM -1

The STREAM-1 series furniture range is a fixed system. Solid worktops are fixed on pedestal units. This system does not require any framework and suits permanent benching particularly with services.

This system is suitable for laboratories that face few and minor changes during their lifetime. Under bench cabinets can be used to support worktops, without need to use C' or H-type frames

Frameless cabinetry is sometimes called "full access" cabinetry because it offers greater accessibility by eliminating the face frame. This makes the cabinets seem larger and may make them more functional and easier to use

**BONLAB STREAM- 1 series are most frequently used system in laboratories..**

Our team's full experience has been applied to our benching systems – engineered to meet the demanding specifications of modern laboratories in all sectors.

